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- Guide sleeve system with matching implant burs suitable for all commonly used implant systems and implant planning programs
- Closed and open molar guide sleeves
- 3D-Laboratory Kit with basic sleeve equipment and template drills for outer sleeves made of titanium
- 3D-Fixation Kit to ensure the best possible fixation of the guide template
- Flapless surgery equipment

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<sup>1</sup> 205= RA L, 206=RA XL <sup>2</sup> Largest working part diameter in 1/10 mm

Hautstanzen  
Biopsy punches  
Punzones de mucosa

Parallel-Pins  
Parallel pins  
Pun guía

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<sup>1</sup> 204=RA <sup>2</sup> Largest working part diameter in 1/10 mm

Standard-Hülsen, geschlossen (je 4 Stück)  
Standard sleeves, closed (4 p. of each size)  
Manguitos estándar, cerrados (4 unid.)

Molare Führungshülsen, offen (je 4 Stück)  
Molar guide sleeves, open (4 p. of each size)  
Manguitos de guía molares, abiertos (4 unid.)

CE

Fig.	HN001	HN002	HN003	HN004	HN005	HN006	HN007	HN008	HN009
Length mm	6,0	6,0	6,0	6,0	6,0	5,5	5,5	5,5	6,0
Internal diameter	2,0	2,5	3,0	3,5	4,0	2,0	3,0	4,0	5,0
External diameter	4,0	4,0	4,0	4,0	4,5	5,0	5,0	5,0	6,0

Internal diameter External diameter

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<sup>1</sup> 205=RA L <sup>2</sup> Largest working part diameter in 1/10 mm

Kugelpfanker (je 10 Stück)  
Ball-head anchors (10 p. of each size)  
Anclajes de bola (10 unid.)

Schraubendreher  
Screw driver  
Destornillador




Spiralbohrer  
Twist drill  
Taladros espirales

CE 0044

Fig.	HN023	HN024	TC05W	UG00W	HN019
Shank <sup>1</sup>			204		205
Size <sup>2</sup>					010
Length mm	8,0	11,0			14,0
Internal diameter	1,0	1,0			
External diameter	1,3	1,3			1,0

<sup>1</sup> 204=RA, 205=RA L <sup>2</sup> Largest working part diameter in 1/10 mm

Internal diameter External diameter

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

85FL069-0409

**Hager & Meisinger GmbH**  
Hansemannstr. 10  
41468 Neuss • Germany  
Tel.: +49 (0) 21 31-20 12 0  
Fax: +49 (0) 21 31-20 12 222  
Internet: www.meisinger.de  
www.bone-management.de  
E-mail: info@meisinger.de

**Meisinger USA, L.L.C.**  
7442 South Tucson Way, Suite 130  
Centennial, Colorado 80112 • USA  
Tel.: +1 (303) 268-5400  
Toll free: +1 (866) 634-7464  
Fax: +1 (303) 268-5407  
Internet: www.meisingerusa.com  
www.bone-management.com  
E-mail: info@meisingerusa.com

MEISINGER  
GERMANY / USA

3D-Navigation-Control  
3D-Navigation-Control Professional

## 3D-Navigation-Control developed with Dr. Nickenig

### 3D-Navigation-Control Professional

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After the use of three-dimensional imaging procedures (e.g. digital volume tomography (DVT), computer tomography (CT)) and after the virtual planning of the implant position using a planning software, precise guide sleeve systems are required to produce the 3D surgical guide templates (using static navigation); during implant surgery, these systems ensure the guiding of the implant cutter in the axis of the virtually planned implant. The 3D-Navigation-Control is a guide sleeve system that is suitable for all commonly used implant systems (and planning programs) and permits the use of surgical guide templates even in tight molar situations. The special tool requirements of flapless surgery are taken into account as much as the (optional) fixing of the surgical guide template via a 3D-Fixation Kit. Equipped with the 3D-Laboratory Kit, the dental technician is able to accurately install the sleeves in the template.

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The closed guide sleeves (standard sleeves) are primarily used for guiding the implant bur in the front molar and premolar area. For the posterior molar area, laterally open molar sleeves that facilitate the lateral placement of the implant bur may be used as an alternative. Depending on the mouth opening and / or the interalveolar distance, it is also possible to use standard sleeves for the posterior molar area. However, in such a case, a prior clinical check using a hand piece and an inserted implant bur is recommended. The set includes standard sleeves with an inner diameter of 2.0 mm (yellow) / 2.5 mm (red) / 3.0 mm (blue) / 3.5 mm (green) and 4.0 mm (violet). The guide sleeve with an inner diameter of 4.0 mm (outer = 4.5 mm) is intended to be used as an outer sleeve, which is firmly integrated into the

template by the dental technician and serves as receptor sleeve for smaller-diameter guide sleeves. As and when required, it is possible to carry out up to 5 guided drillings. The use of a minimum of three sleeve / drill sequences for one implant with a standard diameter (of around 4.0 mm) is recommended. The placement of an implant via guide sleeves is generally not necessary but may be considered if there is insufficient bone quality.

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The implant burs that match the standard sleeves by colour and inscription coding enable the user to prepare the implant bed irrespective of the implant system used. If cylindrical implants are used and depending on the bone quality, the implant bur that is used last needs to be changed to a different diameter. If a conical preparation of the implant bed is required, following the cylindrical preparation using the 3D-Navigation-Control system, the last implant bur recommended for conical implants of the respective implant system should be used last. The implant burs that are included in the 3D-Navigation-Control system are available in three different lengths. The 17.0 mm implant bur is meant to be used as a reserve bur if there is reduced space (only for Ø 2.0 mm). The 20.0 mm implant bur is the standard length that is used most. The 23.0 mm implant bur is primarily intended for single-tooth implants in the front tooth area. Due to the depth markings on the implant bur, the user can control proceedings when the template is removed. After the depth stop has been calculated, the implant bur length (17.0 mm / 20.0 mm / 23.0 mm) indicated in the set equals the exit length of the implant bur from the contra-angle hand piece (use own hand piece to check beforehand!). The bur head length of 0.5 mm must be added. During the virtual planning in the planning program, there should never be less than a minimum distance of 2 to 3 mm from adjacent structures! A "dapping" use of the implant bur as well as the repeated rinsing of the implant bed with a sterile rinsing liquid (e.g. each time the guide sleeves are changed) is recommended. At the latest, the implant bur should be exchanged after it has been used five times.

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The molar guide sleeves are available with inner diameters of 2.0 mm / 3.0 mm / 4.0 mm and with 5.0 mm as an outer sleeve. The lateral opening allows the lateral insertion of the implant bur in cases where the interalveolar distance is greatly reduced or the mouth opening is insufficient. The guiding of the implant bur is reduced in comparison with the standard sleeve.

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The basic sleeve equipment of the 3D-Laboratory Kit includes three standard and three molar sleeves in all diameters. Two template drills for the respective outer sleeve of the standard and / or molar sleeves are enclosed. A screw-down drill extension (BV040 and SW150) allows an adjustment to various scenarios.

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When dealing with edentulous maxilla or very few remaining teeth, it is possible that the fixation of the guiding template during the implant preparation is insufficient. In such cases, it may be necessary to stabilise the template on ball anchors that have temporarily been inserted into the bone. The support provided by ball anchors is based on the patrizie-matrix principle. In the case of a toothless lower jaw, a minimum of three ball anchors is inserted on both sides of the jaw in the retromolar gap and in the (central) anterior tooth area of the patient who has been put under a local anaesthetic. The set comprises a pilot drill as well as a screw driver – to be assembled - which is used to screw in the ball anchor. Depending on the thickness of the mucous membrane and the bone height and/or the distances to adjacent bone structures, a ball anchor with a length of 10.0 mm or 13.0 mm may be used (if necessary as part of a flapless procedure). The insertion of the matrices into the template is generally carried out inside the mouth. Prior to the polymerisation process, a dental dam should be placed on the ball section (patrizie) for protection purposes. The 3D-Fixation Kit is generally used before X-ray images are taken as the position of the template must be the same during the X-ray procedure and surgery. As the anchor is only meant to be in place for a short time, the period between the day when the X-rays are taken and the day when the actual surgery is carried out should be no longer than one week. The ball anchors are not suitable to support the dental prosthesis.

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As a rule, flapless surgery requires prior 3D diagnostics with the planning results being transferred to a guide template. When carrying out flapless surgery, strict compliance with a standard protocol is necessary. A "dapping" use of the implant bur as well as the repeated rinsing of the implant bed

with a sterile rinsing liquid (e.g. each time the guide sleeves are changed) is recommended. The tissue punches that come in different diameters match the different implant diameters. It is recommended to wait with the punching until the first (2.0 mm) drill has been carried out as this will ensure the best punching position. As due to the closed procedure, the direct surgery site cannot be viewed, the use of parallel pins and a measuring gauge to keep a check on and/or protect adjacent bone structures is of great importance.

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##### (Einstellung Tiefenstopp)

##### User-defined adjustment of guide sleeves in a planning program

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If the planning software allows a user-defined (individual) height adjustment of the guide sleeves (depth stop), the following parameters are required for the 3D-Navigation-Control system:

drilling depth:	17.0 mm or 20.0 mm or 23.0 mm
sleeve length:	6.0 mm for standard guide sleeves and
	5.0 mm for molar sleeves.

The length of 0.5 mm for the bur head must always be added (i.e. the total length of the 20.0 mm implant bur: 20.0 mm + 0.5 mm). Hence, it must be pointed out that for safety reasons a minimum distance of 2-3 mm to adjacent structures should be kept during the virtual planning stage in a planning program. All surgery involving the use of 3D-Navigation-Control requires clinical monitoring in accordance with the findings and standards governing implant surgery in particular prior to and during the actual intervention.

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