**Diamantinstrumente**  
Diamond instruments  
Instrumentos de diamante

CE 0044



Fig.	801	801	801	801	HM 1S	HM 1S	HM 1S
Shank <sup>1</sup>	204	204	204	206	204	204	205
Size <sup>2</sup>	023	027	033	021	023	027	023

<sup>1</sup>204=RA, 205=RA L, 206=RA XL <sup>2</sup> Largest working part diameter in 1/10 mm

**Hartmetallbohrer**  
Tungsten carbide burs  
Fresas de carburo

**Handinstrumente**  
Hand instruments  
Instrumentos de mano

Fig.	TL1	TL2	TL3	TL4	TL5	TL6	TL7	TL8	TL9
Name	Sinuslift 1	Sinuslift 2	Sinuslift 3	Sinuslift 4	Sinuslift 5	Sinuslift 6	Membrane Positioning	Condenser	Scissor
Length mm	173.0	173.0	172.0	174.0	175.0	173.0	170.0	181.0	116.0

Art.-No.: BL00

Bone Management® is a registered trademark of Hager & Meisinger GmbH, Germany

**External-Lift-Control**  
Direct Sinuslift System

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
The Bone Management® System External-Lift-Control facilitates easy execution of an external sinus lift in case of insufficient bone height in the lateral maxillary region. To provide adequate bone volume for implant surgery, the interior mucosa of the maxillary sinus (Schneider's membrane) is lifted carefully in the maxillary sinus. This creates additional space for new bone formation which can be filled with autologous bone grafts or bone replacement material.

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In a first step, the gingiva is folded back in the maxillary sinus by appropriate incision to uncover the maxillary bone for further treatment. Using one of the ball-shaped diamonds or hard metal drills included in the set, a window to the sinus is prepared carefully and the bone thinned down to the membrane in this region. Please note that this window is to be folded in an upward direction.

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
Then proceed in the maxillary sinus using one of the angled manual sinus lift instruments. Depending on circumstances, use the different angles and sizes of the manual instruments to carefully separate the mucosa of the maxillary sinus (Schneider's membrane) from the floor of the maxillary sinus without damaging it and thus lifting it upwards.

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
The previously opened bone window in the maxillary sinus wall can now be transferred upwards carefully in the maxillary sinus together with already separated mucosal membrane.

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The created cavity in the maxillary sinus can now be loosely filled with appropriate graft material. This should support the previously lifted bone window. The inserted graft can now be condensed carefully using the Condenser TL8.

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
In case of direct implantation, the implant is placed in the newly created bone site and the remaining open bone access is covered with a membrane, if required, using the TL7 manual instrument „Membrane Positioning“. Then the mucoperiosteal flap is sutured saliva-tight.

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After an appropriate healing period, the implant(s) is/are uncovered followed by appropriate prosthetic restoration. In the two-stage process, implantation is performed only after complete ossification of the graft.

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- Optimal system for the easy execution of an external sinus lift
- Controlled minimal invasive procedure
- Ideal condensation of the implant site

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### General instructions:

- All products delivered are unsterile, therefore, before initial and each further potential application, products should be treated (cleaning / disinfection / sterilization).

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