Fix-Control

Snap-On Anchorage System

Fig.	186RF	E2001	E2003	E2004	E2005	CA0KG	CAFCK	CA00S	I00FC	FC0AK
Shank ¹	204	205	205	205						
Size	018	022	032	034	035					
Length mm	12,0	12,0	12,0	12,0	12,0				12,0	3,5
Minimal diameter			2,24	2,44	2,54				2,54	
External diameter		2,2	3,20	3,40	3,50				3,50	

¹ 204=RA, 205=RA L


Minimal diameter

External diameter

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- Optimum fitting for dentures - simple, cost-effective and can be fitted by every dentist
- Immediate improvement of the denture fit due to softloading
- Implant has optimum initial stability due to high degree of bone compacting, usable for both upper and lower jaws
- No need for complex additional denture components because of the one-piece implant
- Where friction causes loosening, the keyway or O-ring can be easily replaced
- Can be easily combined with telescopic crowns and continuous beams
- Particularly suitable for older patients because the implant just requires normal oral hygiene

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Fix-Control is a one-piece dental implant system made of pure titanium (Grade 4) which can be employed by any dentist to provide a rapid, simple and cost-effective way of enabling patients to retain the same degree of denture comfort even after the last telescopic crown or natural tooth has been lost. Even in cases where the loss of teeth lies some time in the past and it is hardly possible to guarantee that the dentures will stay in position, Fix-Control provides the dentist with the means to improve the patient's quality of life considerably using relatively simple methods and ensure the dentures are firmly held in place again.

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
The starting situation in our example here is a jaw with no teeth left at all into which two dental implants are to be fitted to fix the dentures. After the last telescopic crown, the dentures could not be held in place.

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
In order to drill through the corticalis, which is often extremely hard, it is important to use the initial bur (triangular) first to prevent the pilot bur slipping laterally.

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
Subsequently, the 022 diameter pilot bur will create a preliminary hole slightly under size for orientation and monitoring the angle of the drill (drilling depth up to red colour mark).

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
Use the conical expansion bur I to create a subdimensional cavity. This bur will create a cavity 25% tighter than the implant. If bone density is class D1 or D2, also use conical expansion bur II. This bur creates a subdimensional cavity which is only 5% tighter than the implant. The drill depth depends on the subsequent application depth of the implant. With a gingiva thickness of approx. 2 mm the spreader is used up to the end of the working part (12 mm). Should the shoulder height of the implant lie level with the bone ridge (in gingiva thinner than 2 mm), a drilling with the spreader up to the red colour marking is necessary.

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
The spreader appropriate to the intended implant compacts the bone tissue subsequently and prepares it for the Fix-Control implant. The spreader with the help of

the carrier CA0KG and the driver for carrier is screwed cautiously exactly into the depth position. The implant is not inserted into the cavity until the spreader has taken up the ideal position.

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
In the next stage, the Fix-Control implant is carefully screwed into place using the carrier CAFCK and driver (avoid overwinding). During the healing phase, the existing dentures are drilled out from underneath to form an oversize hole into which the ball head of the implant will easily fit and the dentures are bedded in a plastic masse which remains soft to prevent unnecessary pressure being exerted on the implant.

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inserted into the dentures at right angles to the jaw ridge irrespective of the positioning of the implant.

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General instructions:

- The implants and corresponding keyways can also be ordered separately under the following order no. 100FC. All instruments delivered are unsterile therefore, before initial and each further potential application to be treated (cleaning / disinfection / sterilisation).

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