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- Ridge widening of the mandible by selective spreading and angle modulation of the implant bed
- Controlled, minimally invasive procedures
- Maximal conservation of existing bone substance
- Suitable for preliminary implant cavity preparation for all current implant systems
- Optimal in combination with Split-Control

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
<sup>2</sup> Largest working part diameter in 1/10 mm ↗ Minimal diameter ↘ External diameter

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Horizontal-Control Plus is an angle modulation system that, with its well-proven and unique instruments, provides for simple and minimally invasive ridge widening. This is achieved due to the angle modulation of the implant axis given by the actual bone supply by erecting and displacing of the cortical plate, and at the same time, condensing the bone structure. In this way, the mandible can be optimally prepared for the use of all current implant systems.

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Using the cone cranks – Bone Splitting with angle modulation

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Indications for using Horizontal-Control are "pear-shaped" atrophied or traumatized segments of the mandibular alveolar ridge where a conventional implant is not possible without corresponding bone management.

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First, an undersized pilot hole is made with the depth-marked pilot bur which determines the depth and direction of the implant bed (note diameter and depth setting).

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In the next step, the expansion burs create the bone cavity coordinated with the cone levers that will be used.

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By carefully screwing in a spreader from one of the Bone Management® Split-Control systems, the implant opening is widened in the area of the spongiosa,

gently compressing the adjacent bone tissue in the process.

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Screw one of the horizontal cone cranks – depending on accessibility in the oral cavity – onto the handle and then attach the suitable opening on the spreader. By gradually changing the orientation of the instruments to the angle modulation of the implant axis in the direction of the bone deficit, the alveolar ridge can be widened by up to 8 mm.

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Depending on the indication, this process can be repeated with a spreader with a larger diameter so that the implant bed is properly prepared for inserting the implant.

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The implant inserted immediately following this step holds the expanded gap open and stabilizes the cortical walls.

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If the existing conditions in the jaw do not allow sufficient expansion, an internal vertical osteotomy can be performed in the area of the implant opening from distal and medial upwards (without making a flap). This creates a predetermined

breaking point that breaks open the next time the opening is expanded. In this way, the cortical bone is shifted slightly outwards so the alveolar ridge can be spread to the required width.

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Using one of the two wedge modulators, now carefully bend the moveable cortical plate towards the bone defect.

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