Interim-Control

Bone Auxiliary Implant System

	Vorkörner Initial bur Fresa inicial	Diamantierte Säge Diamond saw Disco de separar	Pilotbohrer Pilot bur Fresa piloto	Parallel-Pin Parallel pin Pun guía	Erweiterer Expansion burs Ensanchadores	Spreader Spreader Spreader	Carrier Carrier Carrier	Schlüssel für Carrier Driver for carrier llave de carrier	Hilfsimplantat Auxiliary implant Implante auxiliar (4 St. / 4 p. / 4 pcs.)		
Fig.	186RF	231DC	A1001	A1002	A1003	A1004	A1005	CA0KG	CA0HI	CA00S	100HI*
Shank ¹	204	204	206	206	206	206	206				
Size ²	018	070	010	009	023	025	027				
Length mm	12,0	0,3	15,0	18,0	15,0	15,0	15,0				
▲					1,03	1,23	1,43				
◻					2,30	2,50	2,70				

¹ 204=RA, 206=RA XL

² Largest working part diameter in 1/10

▲ Minimal diameter

◻ External diameter

* With the reuse of disposable products the risk of infection cannot be excluded and a risk-free functional safety cannot be guaranteed

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- Fixed temporary prosthetic care – easy, economical, and comfortable
- Suitable for all bone resorption classes (1–5) and bone densities DI–D3
- High primary stability through bone-condensing
- Suitable for preliminary implant cavity preparation of all current dental implant systems
- Suitable for definitive care of small incision gaps in the lower jaw

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The auxiliary implants can also be ordered under the following Art.-No. 100HI.

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The Bone Auxiliary Implant System, suitable for all bone resorption classes, provides economical short term care of the jaw with auxiliary implants. These implants facilitate the process of osseointegration of the permanent implants by absorbing their prosthetic stress.

The patient enjoys the advantage of high prosthetic comfort from the very beginning of the post-surgical period. After the phase of osseointegration of the permanent implants, the temporary prosthetic treatment and the auxiliary implants are removed. Should the production of the temporary prosthetic treatment be carried out in a dental laboratory rather than in the mouth itself, the Temporary System also offers special manipulating implants (Art.No.: 39 0000M).

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
In a narrow alveolar ridge situation, the initial bur prevents the pilot bur from slipping. Where necessary, a cortical osteotomy disc will promote bone-spreading.

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
Subsequently, the pilot bur creates a subdimensional pilot drilling for orientation purposes and for control of the angle direction (pay attention to depth of insertion).

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
The parallel pin placed in the osteotomy then checks the direction of drilling both directly and through imaging (x-ray control).

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In the next step, the conical expansion bur I creates a subdimensional cavity which has a volume difference of 25% less in relation to the auxiliary implant. The expansion bur II is only used with bone density D1 and achieves a subdimensional cavity with only a 5% difference to the auxiliary implant.

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
The implementation of the spreader leads to bone condensing and thus prepares the cavity perfectly for subsequent insertion of the auxiliary implant.

The spreader may be screwed in cautiously with help of the carrier CA0KG and, if necessary, the "driver".

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
The auxiliary implant is then inserted gently with help of the carrier CA0H1 and the "driver" (pay attention to depth of insertion). Two (or more) auxiliary implants serve as a basis for a temporary prosthetic treatment. The (permanent) dental implants are thus enhanced in their osseointegration.

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In cases of restricted space such as absence of a lower incisor, the auxiliary implants, due to their tiny dimensions, may also serve as permanent implants.

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General instructions:

- All instruments delivered are unsterile therefore, before initial and each further potential application to be treated (cleaning / disinfection / sterilisation).

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