Diamantbohrer
Diamond bur
Fresa de diamante

Pilotbohrer
Pilot burs
Fresas piloto

Schraubendreher
Screw driver
Atornillador

CE 0044

Fig.	801	TC203	A2001	TC0SW	UG00W
Shank ¹	204	204	204	204	
Size	023	010	013		
Length mm		12,0	12,0		

¹ 204=RA

Schraubendreher
Screw driver
Atornillador
for UG00W

Schrauben (je 10 Stck.)
Screws (10 p. of each size)
Tornillos (10 pzs. c/u)

**ONLY in
Screw Professional**

Schrauben (je 10 Stck.)
Screws (10 p. of each size)
Tornillos (10 pzs. c/u)

Fig.	SWI50	TC0KS*	TC0OS*	TC0LS*	TC20S*
Length mm		7,0	10,0	13,0	10,0
		1,0	1,0	1,0	1,2
		1,3	1,3	1,3	1,5


 Internal diameter

 External diameter

*  With the reuse of disposable products the risk of infection cannot be excluded and a risk-free functional safety cannot be guaranteed

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Art.-No.: BSR00

Art.-No.: BSRPR

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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Screw System Screw System Professional

Bone Fixation

MEISINGER
GERMANY / USA

Screw System Screw System Professional

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The MEISINGER Screw System is an all purpose screw system for easy and safe fixation of bone segments e.g. for horizontal or vertical augmentative grafts. The kit contains screws in two different lengths (7 mm and 10 mm) that are manufactured from special biocompatible titanium alloy. Screw System Professional contains screws in the length 13 mm and with an additional diameter. The Screw System is the optimal supplement to the Bone Management® System Transfer-Control (Art.-No. BTR00).

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To optimise healing, the donor and recipient sites are prepared first with small dotted drillings (blood dots). Then the bone cylinder is placed onto the accordingly prepared bone area and should thereby rest flatly upon the receiving bone. The fixation of the bone cylinder can take place with a hemostat.

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Subsequently, with the thin spiral drill (1.0 mm) a perforation drilling of the bone cylinder and a drilling up to the counter-cortical of the recipient bone is carried out.

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In the next step, only the bone cylinder is perforated with the larger twist drill (1.3 mm).

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If necessary, counter-sink the screw head with the diamond coated round bur.

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
With the help of the inbus key the screw driver is mounted onto the handpiece (see illustration). Then the screw driver leads the bone screw in the suitable length through the drilling in the bone cylinder up to the bone recipient and sensitively fixes it into the receiving bone area. After a healing process of about 3-4 months the screw can be removed.

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- All purpose screw system for simple and safe fixation of bone segments
- Biocompatible screws made from special titanium alloy in two different lengths
- A perfect addition to the Bone Management® Transfer-Control system

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General instructions:

- All products delivered are unsterile therefore, before initial and each further potential application to be treated (cleaning / disinfection / sterilisation).

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