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A further application possibility to the Horizontal Augmentation with the Crest-Control is found in the implementation according to the biological protocol of the Callus-Distraction.

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General instructions:

- The horizontal spreader can be ordered separately under the order number VCD50. All products delivered are unsterile therefore, before initial and each further potential application to be treated (cleaning / disinfection / sterilisation).

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¹ 204=RA, 205=RA L ² Largest working part diameter in 1/10 mm

	Pilotbohrer Pilot burs Fresas piloto		Erweiterer Expansion burs Ensanchadores		Spreader Spreaders Spreader					
Fig.	A1001	D1001	A1003	D1003	A1005	B1005	C1005	D1005	E1005	F1005
Shank ¹	206	206	206	206	206	206	206	206	206	206
Size ²	010	018	023	030	027	029	031	033	035	040
Length mm	15.0	15.0	15.0	15.0	15.0	15.0	15.0	15.0	15.0	15.0
▲					1.43	1.64	1.84	2.05	2.26	2.79
□					2.70	2.90	3.10	3.30	3.50	4.00

¹ 206=RA XL ² Largest working part diameter in 1/10 mm ▲ Minimal diameter □ External diameter

	Ratsche Ratchet Carraca	Eindrehhilfen Jollying aids Soporte para enroscarse			
Fig.	CARA1	CA1RB	CA0RB	CA8RA	CA9RA
Shank ¹				204	205
Size ²		short	long	short	long
Length mm	84.0				

¹ 204=RA, 205=RA L ² Largest working part diameter in 1/10 mm

Horizontal-Spreader (3 Stck.)
Horizontal spreaders (3 p.)
Spreader horizontal (3 pzs.)

Fig.	VCD50
Width	6.0 mm
Length	9.75 mm
Height	10.6 mm
Lifting Height	max. 5 mm

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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
Split-Control is the minimally invasive alternative to osteotomes. Bone Spreading and Bone Condensing with special screw like instruments (spreaders) achieve a controlled and standardised dilatation of horizontally resorbed bone and a gentle densification of cancellous bone. The alveolar ridge can be optimally prepared for the subsequent process of implant insertion while maintaining the existing bone substance without complicated horizontal or vertical bone grafting.

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
To promote Bone Spreading a crestal relief cut with the osteotomy disk (between multi drilling) is implemented.

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
When Bone Spreading the initial bur prevents the pilot bur from slipping. Subsequently the pilot bur produces a subdimensional bone cavity (pay attention to depth of insertion).

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The successive deployment of a series of spreaders with increasing diameters leads to Bone Spreading and condensing (order of diameter and codification of the instruments must be respected). With the help of the ratchet the spreaders may be rotated in cautiously. In bone density of 1–2 the drilling before deploying the spreader is prepared with the corresponding D1, D2 expansion bur.

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
By having widened the implant cavity in this manner it is now possible to place the suitable implant. Furthermore the increased bone rigidity, achieved

by bone condensation, may lead to greater primary stability.

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In the case of Bone Condensing first the pilot bur produces a subdimensional bone cavity (take notice of diameter and depth of the drilling). In a narrow alveolar ridge situation the initial bur prevents the pilot bur from slipping.

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A sub-dimensional bone cavity for bone densities D3/D4 is then created using the conical expansion burs.

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For Bone Condensing purposes of implant cavities of more than 2.7 mm diameter, use the expansion bur A1 and the spreader A1. For implant cavities exceeding 3.3 mm diameter, the expansion bur D1 and the spreader D1 must be used. After placement of the primary spreader (A1 or D1), the following spreaders in ascending order can be used without further drilling with pilot drills. In non-atrophied bone, the use of the spreaders leads to condensation of the bone and horizontally resorbed areas of the alveolar ridge are expanded. With the help of the ratchet, the spreaders may be screwed in cautiously.

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The condensed and spread bone may now provide greater primary stability once an implant has been placed.

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- Minimally invasive alternative to osteotomes
- Controlled dilatation of horizontally resorbed bone
- Gentle densification of cancellous bone
- Suitable for bone resorption classes 1–3
- Provides high primary implant stability through bone condensing
- Suitable for preliminary implant cavity preparation for all current dental implant systems
- Optimal in combination with Horizontal-Control

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Crest-Control is a Horizontal Bone Spreading system which allows for the pre-determined and minimal invasive horizontal widening of the alveolar ridge in particular in the distal lower jaw. With the completely new horizontal spreader the alveolar ridge can be expanded up to a maximum of 5 mm. Subsequently, an implant of all current systems can be used during the mineralisation phase in the widened segment from which the osteointegration due to a highly-potent body growth factor of the existing callus may be accelerated.

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
Indications for the use of the Crest-Control are horizontally resorbed sections e.g. the lower jaw - alveolar ridge, which do not allow the placement of a conventional implant without suitable Bone Management®. In addition, Crest-Control is especially suited when placing several implants side by side.

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
A crestal scoring on the alveolar ridge is carried out first using the diamond spherical instrument. This prevents the following diamond saw from slipping.

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
The diamond saw allows a crestal osteotomy to be created along the scored ridge. This enables the diamond instrument to be used to deepen the crestal osteotomy.

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In the area of the implant opening, with the help of the osteotomy instrument, the crestal cut can be deepened and prepared accordingly for the application of the horizontal spreader.

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