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- Simple, safe transplant of bone cylinders for the restorative construction of a vertically atrophied jawbone using intelligently coordinated instruments
- Prepare a fresh flat contact surface and a central insertion hole for the osteosynthesis screw in a single operation
- Bone cylinder is securely located with no gap
- Optimum site conditions for secure healing
- Considerably reduced operation time due to systematic procedure
- Suitable for preparing the bone for subsequent insertion of all standard forms of intraoral dental implant fixtures
- Optimum combination with Bone Management® set Screw System

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<sup>1</sup>205=RA L    <sup>2</sup>Largest working part diameter in 1/10 mm  
 Internal diameter    External diameter

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
The most important condition for a safe healing of an osseous graft is a congruent  
and fresh contact surface of the implant area. The Transfer-Ring-Control system  
allows this precondition for the healing of vertical plastic supports in a simple and  
controlled manner – while at the same time reducing operation time considerably.  
The Bone Management® set Screw System (order no. CB0SK) is the ideal comple-  
ment to the Transfer-Ring-Control.

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
The initial phase requires the donor region to be  
marked. The hole thus created is used later to en-  
sure perfect centring of the osteosynthesis screw in the implant site, at the same  
time, the surface of the intended bone cylinder is ground flat thus creating a radial  
groove.

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
Next, the trephine is used to cut a bone cylinder of the  
correct diameter from the donor site. Depending on  
the vertical requirements at the site, this cylinder can  
have a height of 4 mm to 7 mm.

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The bone cylinder is carefully removed from the do-  
nor site which can then be filled with collagen.

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If the surface of the cylinder is still uneven, it can be ground flat in a simple and  
controlled manner using the abrasive drill. Similarly, the central guide hole for  
the osteosynthesis screw can be drilled.

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Next, the abrasive bur is used to form a standardized bone site with a flat, fresh  
contact area. The bone cylinder can then be inserted perfectly into the hole.  
Since the instruments used are all designed to complement each other, the  
transplant will be a perfect, gap-free fit.

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Once the bone cylinder has been inserted in the  
site, it is secured using the osteosynthesis screw  
(can be ideally combined with the screws in the  
Bone Management® set Screw System).

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If required, splinters on the upper surface of the bone cylinder can be removed  
using the diamond grinder.

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## General instructions:

- All instruments delivered are unsterile therefore, before initial and each further potential  
application to be treated (cleaning / disinfection / sterilisation).

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