Transfer-Control

Bone Replacing System

CE 0044

Fig.	T229L	T229L	T229L	TC084	TC084	TC084	TC21X	TC21X	TC21X
Shank ¹	104 / 205	104 / 205	104 / 205	204	204	204	104 / 205	104 / 205	104/205
Size ²	040	050	060	004	005	006	040	050	060
Length mm	14,0	14,0	14,0	4,0	5,0	6,0	14,0	14,0	14,0
	4,0	5,0	6,0						
	5,0	6,0	7,0	10,0	10,0	10,0			

¹ 104=HP, 204=RA, 205=RA L

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- Safe and easy transplantation of autogenous bone cylinders with a set of precise tools
- Optimally suitable for horizontal and vertical augmentative grafts
- Predictable application utilizing precise depth markings
- Suitable for preliminary implant cavity preparation for all current dental implant systems

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¹ 204=RA, 205=RA L

²Largest working part diameter in l/10 mm

Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

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
This Bone Replacing System, suitable for use with all current dental implant systems, permits precise and standardised transplantation of autogenous bone cylinders. Different types of bone burs are available according to three groups of diameters. Within each group the instruments are adapted to each other so that the outer diameter of the ablative bur or the diameter of the wheel cutter correlates with the inner diameter of the trephine. A press fit is often produced between the cylindrical bone autograft and the prepared recipient site, which can be strengthened with a fixation screw (sold separately). Such precise fitting of autografts lead to accelerated bone revitalisation and wound-healing, thus creating a perfect site for dental implantation within 3 to 4 months.

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In cases of severe horizontal alveolar ridge resorption (where even bone-spreading techniques cannot be deployed) **horizontal augmentative grafts** are required for compensation purposes.

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Initially an ablative bur (ensure diameter) is chosen and utilized to drill a standardised cavity against the horizontally resorbed alveolar ridge according to the diameter of the subsequently inserted implant. The surgical technique and the properties of the ablative bone bur lead to a hyperaemic environment, which accelerates the wound healing in post transplantation.

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With the appropriate trephine, corresponding to the previously used ablative bur, a bone cylinder congruent to the diameter of the recipient region is extracted e.g. from the retro-molar or other selected donor site. If necessary, the osteotomy disk or diamond can be used to create a horizontal cut allowing the bone cylinder to be taken from the donor site.

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The donor site may be filled with collagen or other grafting material. The intervention is terminated by a current gingival suture.

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The extracted bone cylinder is placed into the recipient bone site to compensate the horizontal bone deficit. A satisfying press fit may render additional fixation unnecessary.

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If a certain mobility of the autograft in the cavity remains, a fixation screw from the Screw System (Art.-No. BSR00, BSRPR) which optimally corresponds with Transfer-Control, may be implemented.

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In cases of severe vertical alveolar ridge resorption, **vertical augmentative grafts** are required for compensation. In this case the appropriate wheel-cutter creates a congruent cavity in the crestal alveolar ridge. One should take note of the correlating trephines to be used. Prepare an autogenous bone graft in the same manner as described for horizontal defects.

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General instructions:

- All products delivered are unsterile therefore, before initial and each further potential application to be treated (cleaning / disinfection / sterilisation).

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