# Master-Pin-Control

## Master-Pin-Control Basic

### Master-Pin-Control Plus

85FLBM08 - 0923

## Pin System for Membrane Fixation

Art.-No.: BMP00, BMPBA, BMPPL

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## Sterilization of the pins

The instruments and pins are supplied non-sterile and need to be cleaned, disinfected and sterilized prior to use.

For proper processing of the articles, please note the instructions for processing (cleaning, disinfection and sterilizing) of medical devices produced by Hager & Meisinger GmbH.

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## Disposal

When disposing of the instruments (at the end of their service life or after the specified shelf life has expired), ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national regulations (e.g. Dual Waste System).

85FLBM08 - 0923

# Application and safety instructions for the MEISINGER Master-Pin-Control

## Master-Pin-Control Basic

### Master-Pin-Control Plus

85FLBM08 - 0923

## Pin System for Membrane Fixation

Art.-No.: BMP00, BMPBA, BMPPL

Master-Pin-Control is designed especially for the fixation of resorbable, non-resorbable and titanium reinforced membranes. The included pins consist of titanium alloy Grade 5 and have a groove which increases the surface of the pins. With their sharp tip and very stable shank, they can easily and precisely be inserted into dense cortical bone.

## Indications for Use

The Bone Management® Master-Pin-Control systems are indicated for the fixation of apical membrane ends at the local bone in order to avoid micromobility of the membrane. They allow for the fixation of resorbable, non-resorbable and titanium-reinforced membranes.

## Utilization

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed.

## Placing of bleeding points

For freshening the bone, bleeding points can be placed in the cortical bone in the region to be augmented using the included twist drills. If necessary, the bone can be pre-drilled with the aid of the initial drills. In case of hard bone, pre-drilling is recommended.

## Application of the augmentation material

According to the respective manufacturer's instructions for use, the augmentation material is applied into the present defect and the selected membrane is adapted to the defect.

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The pins are located in a sterilizable box, which can be opened by a sliding mechanism on the side. With the aid of the pin holder, the pins are taken from the box. For this, the pin holder is pressed vertically from above onto the pin heads. It is important to pay attention to a secure hold of the pins in the pin holder. In order to align the pins straight in the pin holder, the pin holder is optimally held between two fingers when pressing on the pins and then briefly moved counter-clockwise and clockwise once on the pin head in an axial rotary movement. This centres the pins in the pin holder. The pins are not fixed in the box. Therefore, when opening the box as well as when picking up the pins with the pin holder, ensure vibration-free working.

Since the pins are implantable products, the pins should only be picked up using the pin holder to avoid touching them by hand.

## Fixation of the membrane

To fix the membrane over the defect, the pins clamped in the pin holder are applied through the membrane into the cortical bone using a surgical hammer (e.g. HI070) while making sure that aspiration of the pins does not occur. By laterally tilting the pin holder, the holder can be released from the applied pin without negatively affecting the fixation of the pin in the bone. Vertical removal of the pin holder should be avoided as this may cause excessive pulling forces on the pin. In case of very soft bone, the use of the fixing holder is optionally indicated in order to fix the pins in the bone when releasing the pin holder. For this purpose, insert the uncoated working part of the fixing holder between the lamellae of the pin holder while disengaging the pin holder from the pin.

Important: Compression of the membrane must be avoided.

## Removal of the pins

After a clinically confirmed, successful healing period, a small incision of the gingiva is made above the pin to be removed. The pin can then be carefully loosened, for example by using a flat raspatorium, and removed with forceps. With the MP10 pins, the MP1I screwdriver can best be used for loosening. In this case, the pin is loosened by turning it counter-clockwise by using the screwdriver. The pin can then be removed with forceps.

NOTE: When using the screwdriver MP1I, please note that no clamping effect is provided between the screwdriver and the pin MP10. The pins are a disposable product and are therefore not to be used again. The membrane should be handled according to the manufacturer's instructions for use.

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