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Application and safety instructions for the MEISINGER Tenting and Umbrella Screws

Art.-No.: TS035, US050, US070

BSFLBMB32 - 1023

These user instructions apply to all the products listed in the article overview Bone screws. All screws described are made of titanium alloy Grade 5-ELI (titanium grade 23) and have a self-tapping thread. The screw head includes an interface for receiving the insertion tool (hexagon socket).

Indications for Use

The Tenting and Umbrella Bone Screws are designed for the fixation and stabilization of appropriate bone grafts and bone substitutes, as well as for the support of barrier membranes. They are used to restore anatomical structures in dental surgery (osteosynthesis, correction of degenerative illnesses), but are not suitable as substitutes for normal body structures.

Utilization

The bone screws can be used for a variety of surgical techniques depending on the case situation. The choice of technique as well as the compilation of the instrumentation to be used must be appropriately selected and prepared for deployment prior to the start of surgery. Due to the self-tapping thread design, no pilot bore is necessary for the insertion of the screws.

Compatibility

Only Meisinger screwdrivers with hexagon socket (e.g. 31053, TL0T1) may be used for the Tenting and Umbrella Screws. The use of third-party accessory products may impair the function of the screws and excludes any guarantee or replacement.

Postoperative procedure

After clinically and radiographically confirmed complete bone healing and functional restoration, the screws are either removed or left in situ, depending on the case situation and surgical goal. If infection or weakening of the bone tissue occurs as a result of the inserted bone screw the screw should be removed. Explantation is also essential in the case of screw loosening and always in children to avoid growth retardation.

Sterilisation of the screws

The screws are supplied non-sterile and must be prepared before use (cleaning / disinfection / sterilisation). For proper reprocessing of the items, please observe the comments on reprocessing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH.

Disposal

When disposing of the instruments (at the end of their service life or after the specified shelf life has expired), ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national regulations (e.g. Dual Waste System).

Contraindications


Basically, general medical as well as local, absolute and relative contraindications for dental surgical measures must be observed.

Absolute contraindications

- Active infections in or near the area to be augmented as well as local or systemic pathological processes (e.g., symptoms such as fever, local inflammation, abscesses)
- Inadequate bone supply (quality / quantity)
- Illnesses or medication that impair bone metabolism
- Inadequate soft tissue supply, which does not allow complete coverage of the augmented area

Relative contraindications

- Dentoalveolar growth not completed (exception: cases where no dentoalveolar growth is expected, e.g., ectodermal dysplasia) cases where no dentoalveolar growth is expected, e.g., ectodermal dysplasia)
- Known allergies and / or foreign body reaction to the implant material
- Systemic and/or metabolic illnesses or medical treatments that lead to progressive bone deterioration (e.g., cortisone, immunosuppressants, bisphosphonates and diabetes)
- Drug and alcohol abuse
- Lack of patient compliance
- Poor circulation
- Heavy work or active sports
- Mental state which may lead to non-compliance with the medical order
- Highly atrophied jaw

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L₄ (mm) 7,75 8,75 9,75 7,0 9,0 11,0 7,0 9,0 11,0

CE 0044 /

- Nicotine consumption, which may jeopardize the success of the procedure by delaying bone wound healing

IMPORTANT:
Care must be taken to protect the anatomical structures (safety distance min. 2 mm) as well as the course of the adjacent teeth / tooth roots (hazard of damage, infection / resorption). To avoid the development of infections and dehiscence, care should be taken to prevent contact between the graft and adjacent teeth. Cases in which the mechanical stability and mechanical rest of the graft cannot be achieved may compromise successful bone healing.

SEPERATE INSTRUCTIONS FOR SCREWS:
When selecting the screw length as well as the screw placement, the evaluation of the available bone supply by means of X-ray diagnostics must be taken into account. To prophylactically counteract soft tissue perforation or dehiscence, tension-free and complete wound closure must be achieved. The entire augmentation area, including the bone screws, must not be loaded during the healing process. Failure to do so may result in bending, loosening, loss and/or breakage of the screws and this means a negative impact on the augmentation result. Also, excessive physical activity as well as trauma affecting the augmentation site can lead to premature failure of the screws by loosening or breaking. The screws must only be used in accordance with the instructions for use and safety. The complete sinking of the screw can be limited by very hard bone quality. Any kind of modification in the procedure or on the product can lead to impairment of the treatment success or damage to the product. The screws are a single-use product and must therefore not be modified and reused after use. Even screws that have already been used but are visually undamaged may suffer microscopic damage from previous mechanical stress leading to fatigue and product failure. The screws have not been evaluated for safety and compatibility in the MR environment. They were not tested for heating, migration, or image artefacts in the MR environment. The safety of the screws is uncertain in the MR environment. Scanning a patient with this medical device could result in injury to the patient.

RISK OF ASPIRATION: The user must ensure that the screws are secured against aspiration both during the attachment and during its removal.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

General instructions
Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Meisinger seit 1888 CE / CE 0044 Conformity mark Do not use if the packaging is damaged Manufacturing date

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