
Inhalt Transfer-Ring-Control Professional

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


- Simple, safe transplant of bone cylinders for the restorative construction of a vertically atrophied jawbone using intelligently coordinated instruments
- Flexibility due to a great selection of instrument sizes
- Prepare a fresh flat contact surface and a central insertion hole for the osteosynthesis screw in a single operation
- Safe gap-free fit of the bone cylinder offering an ideal condition for a safe wound healing period
- Considerably reduced operation time due to systematic procedure
- Suitable for preparing the bone for subsequent insertion of all standard forms of intraoral dental implant fixtures

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Bone Management® is a registered trademark of the Hager & Meisinger GmbH, Germany

Transfer-Ring-Control Basic Transfer-Ring-Control Professional

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Meisinger since 1888

Transfer-Ring-Control Basic

Transfer-Ring-Control Professional

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The MEISINGER Transfer-Ring-Control systems are indicated for intraoral bone cylinder transplantation in the case of vertical bone defects. With the aid of these systems, the most important conditions for the healing of vertical augmentative grafts, as well as a congruent and fresh bone recipient site can be created in an easy and controlled manner. This system offers a variety of various sizes of instruments which provides high flexibility when removing bone cylinders. Depending on what is needed, various large initial bur trephines, trephines and end milling cutters are available which allow a single-stage as well as simultaneous work. With the help of the additional diamond and tungsten carbide bur, bone splinters can be abraded to precision.

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

The initial phase requires the donor region to be marked with the initial bur trephine. The pilot drilling thus created is used

later to ensure perfect centering of the osteosynthesis screw in the implant site. At the same time, the surface of the intended bone cylinder is ground flat thus creating a radial groove.

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
Next, the trephine is used to cut a bone cylinder of the correct diameter from the donor site. Depending on the

vertical requirements at the site, this cylinder can have a height of 6 mm to 14 mm. The diameter of the bone cylinder can range from 4 mm up to 8 mm. **Trepanning speed may not exceed 300 min⁻¹.**

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With the help of the hand instrument, the bone cylinder is carefully taken from the donor region, and then subsequently the wound area is treated if necessary.

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By using the end milling cutter, surface roughness of the bone cylinder can be milled flat in a controlled manner accordingly. If a fixation of the bone cylinder is necessary, the next step is to drill a central guide hole for the osteosynthesis screw.

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Next, the end milling cutter is used to form a standard sized bone site with a flat, fresh contact area into which

the bone cylinder can be inserted perfectly. Since the instruments used are all designed to complement each other, the transplant will be a perfect, gap-free fit.

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After inserting the bone cylinder into the recipient bone site, this is attached using an appropriate osteosynthesis screw from the Screw System TX (Art.-No. BTX00, BTXPPR), which optimally complements Transfer-Ring-Control.

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If required, splinters on the upper surface of the bone cylinder can be removed using the diamond polisher or the allport bur which are included in the set.

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