



Zweistufig diamantierte Instrumente speziell zur Bearbeitung von silikat- und leucitverstärkten Keramiken.  
 Two-stage diamond coated instruments especially for the treatment of silicate and leucit strengthened ceramics.  
 Instrumento diamantado de dos capas especialmente para trabajos sobre cerámicas reforzadas en silicato y leucita.

rot / orange (mittel), Glätten und Vorpulieren • red / orange (medium), abrading and pre-polishing • rojo/anaranjado (medio), Alisar y prepulido

## DPO01



| Fig.  | Shank | ISO-No.         |      | 1   |
|-------|-------|-----------------|------|-----|
|       |       |                 | L mm | 2,0 |
| DPO01 | HP    | 803 104 303 522 |      | 200 |



## DPO02



| Fig.  | Shank | ISO-No.         |      | 2    |
|-------|-------|-----------------|------|------|
|       |       |                 | L mm | 13,0 |
| DPO02 | HP    | 803 104 243 522 |      | 040  |



## DPO03



| Fig.  | Shank | ISO-No.         |      | 2   |
|-------|-------|-----------------|------|-----|
|       |       |                 | L mm | 2,5 |
| DPO03 | HP    | 803 104 372 522 |      | 110 |



- für alle Silikatkeramiken
- einfache Handhabung
- lange Lebensdauer
- for all silicate ceramics
- easy use
- long lasting
- para todas las cerámicas de silicato
- fácil de utilizar
- vida larga

grau / orange (fein), Hochglanzpolitur • gray / orange (fine), bright lustre polish • gris/anaranjado (final): Pulido brillante

## DPO04



| Fig.  | Shank | ISO-No.         |      | 1   |
|-------|-------|-----------------|------|-----|
|       |       |                 | L mm | 2,0 |
| DPO04 | HP    | 803 104 303 512 |      | 200 |



## DPO05



| Fig.  | Shank | ISO-No.         |      | 1    |
|-------|-------|-----------------|------|------|
|       |       |                 | L mm | 13,0 |
| DPO05 | HP    | 803 104 243 512 |      | 040  |



## DPO06



| Fig.  | Shank | ISO-No.         |      | 1   |
|-------|-------|-----------------|------|-----|
|       |       |                 | L mm | 2,5 |
| DPO06 | HP    | 803 104 372 512 |      | 110 |



Zweistufig diamantierte Instrumente speziell zur Bearbeitung Keramik, CAD/CAM Restaurationen und Oxydkeramiken.  
 Two-stage diamond coated instruments especially for the treatment of ceramics, restorations CAD/CAM and oxide ceramics.  
 Instrumento diamantado de dos capas especialmente indicado para trabajos en cerámica, restauraciones en CAD/CAM y cerámicas d'óxido.



hellblau / orange (mittel), Glätten und Vorpolieren • light blue / orange (medium), abrading and pre-polishing • Azul claro/anaranjado (medio): Alisar y prepulido

## DCA01



| Fig.  | Shank | ISO-No.                |      | 1   |
|-------|-------|------------------------|------|-----|
|       |       |                        | L mm | 2,0 |
| DCA01 | HP    | 803 <b>104</b> 304 524 |      | 260 |



## DCA02



| Fig.  | Shank | ISO-No.                |      | 1    |
|-------|-------|------------------------|------|------|
|       |       |                        | L mm | 14,0 |
| DCA02 | HP    | 803 <b>104</b> 292 524 |      | 073  |



## DCA03



| Fig.  | Shank | ISO-No.                |      | 1   |
|-------|-------|------------------------|------|-----|
|       |       |                        | L mm | 2,0 |
| DCA03 | HP    | 803 <b>104</b> 303 524 |      | 200 |



## DCA04



| Fig.  | Shank | ISO-No.                |      | 2    |
|-------|-------|------------------------|------|------|
|       |       |                        | L mm | 13,0 |
| DCA04 | HP    | 803 <b>104</b> 243 524 |      | 040  |




## DCA05



| Fig.  | Shank | ISO-No.                |      | 2   |
|-------|-------|------------------------|------|-----|
|       |       |                        | L mm | 2,5 |
| DCA05 | HP    | 803 <b>104</b> 372 524 |      | 110 |



## DCA06



| Fig.  | Shank | ISO-No.                |      | 1   |
|-------|-------|------------------------|------|-----|
|       |       |                        | L mm | 2,5 |
| DCA06 | HP    | 803 <b>900</b> 373 524 |      | 170 |



- speziell für Oxydkeramiken
- zum Anpassen von CAD/CAM Restaurationen
- als Korrekturwerkzeug bei Glanzbrandfehlern

- especially for oxide ceramics
- for the fitting of restorations CAD/CAM as an error
- correction instrument on glaze firing defects

- especialmente para las cerámicas d'óxido
- correcciones de restauraciones CAD/CAM
- herramienta de corrección para los errores de quemado final

braun / orange (fein), Hochglanzpolitur • brown / orange (fine), bright lustre polish • Marrón/anaranjado (fino): Pulido brillante

## DCA07



| Fig.  | Shank | ISO-No.                |      | 1   |
|-------|-------|------------------------|------|-----|
|       |       |                        | L mm | 2,0 |
| DCA07 | HP    | 803 <b>104</b> 304 514 |      | 260 |



## DCA08



| Fig.  | Shank | ISO-No.                |      | 1    |
|-------|-------|------------------------|------|------|
|       |       |                        | L mm | 14,0 |
| DCA08 | HP    | 803 <b>104</b> 292 514 |      | 073  |



## DCA09



| Fig.  | Shank | ISO-No.                |      | 1   |
|-------|-------|------------------------|------|-----|
|       |       |                        | L mm | 2,0 |
| DCA09 | HP    | 803 <b>104</b> 303 514 |      | 200 |



## DCA10



| Fig.  | Shank | ISO-No.                |      | 2    |
|-------|-------|------------------------|------|------|
|       |       |                        | L mm | 13,0 |
| DCA10 | HP    | 803 <b>104</b> 243 514 |      | 040  |




## DCA11



| Fig.  | Shank | ISO-No.                |      | 2   |
|-------|-------|------------------------|------|-----|
|       |       |                        | L mm | 2,5 |
| DCA11 | HP    | 803 <b>104</b> 372 514 |      | 110 |



## DCA12



| Fig.  | Shank | ISO-No.                |      | 1   |
|-------|-------|------------------------|------|-----|
|       |       |                        | L mm | 2,5 |
| DCA12 | HP    | 803 <b>900</b> 373 514 |      | 170 |

