

# Prophylaxe



**MEISINGER**  
GERMANY / USA


## Prophy-Cups / Single-use system

Polishing even inside the gingival pocket, easy and cost efficient use.



## Prophy-Clean

Scale Removal, precise removal of tooth discolorations, easy cleaning of areas difficult to reach, no damages to enamel.



Art. No.	PHA02	PS003	PSS04	PCL01
Packaging size	100 pcs	100 pcs	100 pcs	5 pcs
Color	blue	green	yellow	white
Hardness	hard	soft	super soft	hard
Medical device class	Ila	Ila	Ila	Ila
Size mm	7 x 10	7 x 10	7 x 10	2,35 x 9
ISO-No.:	653 204 034 514 070	653 204 034 512 070	653 204 034 511 070	633 204 172 536 023
RPM	10.000 min <sup>-1</sup>	10.000 min <sup>-1</sup>	10.000 min <sup>-1</sup>	10.000 min <sup>-1</sup>

## Prophy-Stripes

The stripes are made of high flexible polyester and are one-side coated with aluminium oxide particles.



Art. No.	PST01	PST02
Packaging size	50 pcs	50 pcs
Hardness / color	coarse / blue - medium / red	fine / yellow - extra fine / white
Thickness mm	0,15	0,15
Size mm	4 x 150	4 x 150



### Hager & Meisinger GmbH

Hansemanstr. 10  
41468 Neuss • Germany  
Tel.: +49 (0) 21 31-20 120  
Fax: +49 (0) 21 31-20 12222  
Internet: [www.meisinger.de](http://www.meisinger.de)  
[www.bone-management.de](http://www.bone-management.de)  
E-mail: [info@meisinger.de](mailto:info@meisinger.de)

