

new products

ACRYLIC CUTTERS

Hager & Meisinger offers the Acrylic Master Kit with new facet cutters. Described as exceptionally resistant, the reinforced, hard-metal instruments have wide cutting edges and ample chip spacing for maximum expulsion of material. The drills from this set are particularly suited for all permanent soft acrylic resins and other plastics, as well as for gyp-

sum. They remove material quickly and evenly, which makes subsequent polishing easier.

Hager & Meisinger GmbH, Germany
(49) 21312012160
www.meisinger.de
Select 12



DPRE Ausgabe März 2006