Diamond Disc Assortment

The right discs for every application

MEISINGER diamond discs are suitable for precise cutting and separation. Whether plastic, ceramics, zirconia or plaster - the comprehensive assortment contains the right discs for every material.
Together they all have a tendency to perfection: a precise diamond coating, very long service life and vibration-free work characterize the diamond discs. Almost all instruments are available in different diameters.



Made in Germany


At a glance

- Safe application
- Precise diamond coating
- Long service life
- Vibration-free working

Our sales force would be pleased to answer your questions: +49 2131 2012-303

Meisinger since 1888

943DC extra fine






Fig.	Shank	Ref.-No.	1	1
		L mm	0,15	0,15
943DC	HP	806 **104** 361 504	080	100

365D fine

New







Fig.	Shank	Ref.-No.	1	1
		L mm	0,3	0,35
365D	HP	806 **104** 333 514	300	450*

910D medium






	Shank	Ref.-No.	1	1
		L mm	0,5	0,5
910D	HP	806 **104** 321 524	190	220
	unmounted	806 **900** 321 524	190	220

911D medium

911DF fine

New






Fig.	Shank	Ref.-No.	1
		L mm	0,5
911D	HP	806 **104** 332 524	220*
	unmounted	806 **900** 332 524	220*
		L mm	0,45
911DF	HP	806 **104** 332 514	220
	unmounted	806 **900** 332 514	220

915D medium

915DF fine
915DC extra fine

New






Fig.	Shank	Ref.-No.	1	1
		L mm	0,25	0,25
915D	HP	806 **104** 355 524	190*	220*
	unmounted	806 **900** 355 524	190*	220*
		L mm	0,2	0,2
915DF	HP	806 **104** 355 514	190	220
		L mm	0,15	0,15
915DC	HP	806 **104** 355 504	190	220

916D medium

916DH super coarse
916DG coarse
916DF fine






Fig.	Shank	Ref.-No.	1	1
		L mm	0,5	0,5
916D	HP	806 **104** 327 524	190	220
	unmounted	806 **900** 327 524	190	220

918D medium






Fig.	Shank	Ref.-No.	1
		L mm	0,5
918D	HP	806 **104** 335 524	220
	unmounted	806 **900** 335 524	220

921DF fine

921DC extra fine

New






Fig.	Shank	Ref.-No.	1	1
		L mm		0,2
921DF	HP	806 **104** 358 514		220*
	unmounted	806 **900** 358 514		220*
		L mm	0,15	0,15
921DC	HP	806 **104** 358 504	190	220
	unmounted	806 **900** 358 504	190	220

930D medium

930DF fine






Fig.	Shank	Ref.-No.	1	1
		L mm	0,3	0,3
930D	HP	806 **104** 345 524		220
	unmounted	806 **900** 345 524		220
		L mm	0,25	0,25
930DF	HP	806 **104** 345 514	190	220
	unmounted	806 **900** 345 514	190	220

932D medium

932DF fine






Fig.	Shank	Ref.-No.	1	1
		L mm	0,3	0,3
932D	HP	806 **104** 350 524		220
	unmounted	806 **900** 350 524		220
		L mm	0,25	0,25
932DF	HP	806 **104** 350 514	190	220
	unmounted	806 **900** 350 514	190	220

935D medium

935DF fine

New






Fig.	Shank	Ref.-No.	1	1
		L mm	0,3	0,3
935D	HP	806 **104** 340 524		220*
	unmounted	806 **900** 340 524		220*
935DF	HP	806 **104** 340 514	190	220
	unmounted	806 **900** 340 514	190	220

947DF fine

947DC extra fine

New






Fig.	Shank	Ref.-No.	1	1
		L mm	0,2	0,2
947DF	HP	806 **104** 389 514	190	220
	unmounted	806 **900** 389 514	190	220
		L mm		0,15
947DC	HP	806 **104** 389 504		220*
	unmounted	806 **900** 389 504		220*

940F fine






Fig.	Shank	Ref.-No.	1	1	1
		L mm	0,2	0,2	0,2
940F	HP	806 **104** 377 514	180	200	220
	unmounted	806 **900** 377 514	180	200	

401DF fine

New






Fig.	Shank	Ref.-No.	1
		L mm	0,15
401DF	HP	806 **104** 400 514	220
	unmounted	806 **900** 400 514	220

Recommended use: Diamond Discs

Diamond Disc				Remove plastic	Seperate plastic	Remove ceramic	Seperate ceramic	Zirconia	Plaster
Figure	L	Ø	Handling						
910D	0,5	190 220	rigid			✓		✓	
911D	0,5	220	rigid	✓		✓		✓	
	0,45	220						✓	
915D	0,25	190 220	super flexible		✓		✓		
	0,2	190 220							
	0,15	190 220							
916D	0,5	190 220	rigid			✓			
918D	0,5	220	rigid	✓		✓			
921D	0,2	220	super flexible		✓		✓		
	0,15	190 220							
930D	0,3	220	flexible			✓		✓	
	0,25	190 220			✓		✓		
932D	0,3	220	flexible	✓		✓	✓	✓	
	0,25	190 220			✓			✓	
935D	0,3	220	flexible	✓	✓	✓	✓	✓	
	0,3	190 220							
940F	0,2	180 200 220	super flexible		✓		✓		✓
943D	0,15	080 100	super flexible		✓		✓		
947D	0,2	190 220	super flexible		✓		✓		
	0,15	220							
401D	0,15	220	flexible		✓		✓		
365D	0,3	300	flexible						✓
	0,35	450							

✓ = particularly suitable ✓ = suitable

Meisinger USA. L.L.C.
10150 E. Easter Avenue • Centennial, Colorado 80112 • USA
Phone: +1 (303) 268-5400 • Fax: +1 (303) 268-5407
Toll free: +1 (866) 634-7464
info@meisingerusa.com • www.meisingerusa.com

Please observe the general instructions for the use and safety of MEISINGER products in the medical field as well as the notes on reprocessing (cleaning, disinfection and sterilization) Hager & Meisinger GmbH medical devices.

85FL116EN - 0317